

## Power HotStart *Taq* DNA Polymerase

**Cat. No:** PR1601    200 U (40 µl)  
                  PR1602    1000 U (200 µl)

**Storage:** Up to one year at -20°C

### **Description:**

Hot Start *Taq* DNA Polymerase is a recombinant *Taq* DNA polymerase, which has been chemically modified by the addition of heat-labile blocking groups to its amino acid residues. The enzyme is inactive at room temperature, avoiding extension of non-specifically annealed primers or primer dimers and providing higher specificity of DNA amplification. The functional activity of the enzyme is recovery during short 4-minute incubation at 95°C. The activated enzyme possesses the same functionality as *Taq* DNA polymerase: catalyzes 5'→3' synthesis of DNA, having no detectable 3'→5' proofreading exonuclease activity, but possesses low 5'→3' exonuclease activity. It exhibits deoxynucleotidyl transferase activity, which frequently results in the addition of extra adenines at the 3'-end of PCR products. Before activation, the two activities are not detectable.

### **Unit definition :**

One unit is defined as the amount of enzyme required to incorporate 10 nm isotope labeling dNTP into insoluble acid material at 72°C in 30 minutes

### Reaction Setup

Hotstart <i>Taq</i> DNA Polymerase (5 U/μl)	0.5-1 μl
10×PCR Buffer (Mg <sup>2+</sup> plus)	5 μl
dNTP Mixture(2.5 mM each)	4 μl
Template DNA	< 1 g
Primer 1 (20 μM)	1 μl
Primer 2 (20 μM)	1 μl
Sterile Water	Up to 50 μl

### Thermal cycler condition

Cycle Numbers	Step	Temperature	Time
1	1	95°C	2-5 min
30-35	1	95°C	30 s
	2	50-60°C	30 s
	3	72°C	1 min/1-2kb
1	1	72°C	10 min

### Notes

- Please optimize the reagent concentrations, conditions, and parameters according to the experiment.
- Make PCR reaction mixture on ice to avoid non-specific reaction.