

*Note: for laboratory research use only*

## **Whole Blood DNA Extraction Kit (Spin-column)**

**Cat. #: DP1801 (50preps)  
DP1802 (100preps)**



**Signalway Biotechnology**

---

# **I. Kit Content, Storage, and Stability:**

| <b>Component</b>                         | <b>Storage</b> | <b>50 preps<br/>(DP1801)</b>            | <b>100 preps<br/>(DP1802)</b>            |
|------------------------------------------|----------------|-----------------------------------------|------------------------------------------|
| Buffer BB                                | RT             | 15 ml                                   | 30 ml                                    |
| Binding Buffer CB                        | RT             | 15 ml                                   | 30 ml                                    |
| Inhibitor<br><br>Removing Buffer<br>(IR) | RT             | 27 ml                                   | 50 ml                                    |
| Washing Buffer<br><br>(WB)               | RT             | 15 ml                                   | 25 ml                                    |
|                                          |                | <i>Add 60 ml ethanol<br/>before use</i> | <i>Add 100 ml ethanol<br/>before use</i> |
| Eluting Buffer EB                        | RT             | 15 ml                                   | 20 ml                                    |
| Isopropanol                              | RT             | 7 ml                                    | 15 ml                                    |
| Protease K (20<br>mg/ml)                 | -20°C          | 20 mg<br><br>Dry powder                 | 20 mg×2<br><br>Dry powder                |
| Spin-column AC                           | RT             | 50 pcs                                  | 100 pcs                                  |
| Collection Tube<br>(2ml)                 | RT             | 50 pcs                                  | 100 pcs                                  |

*All reagents are stable for 12 months at RT if stored properly.*

## II. Notes :

1. *Add ethanol to Buffer WB before use, mix adequately, and then check the box on the label showing ethanol was added!*
2. Precipitate may occur in buffers under low temperature, and it will affect DNA yield. Incubate to 37°C for a few minutes **until buffer is clear**, and then let the buffer cool to RT before use.
3. *Proteinase K is provided in freeze-dried powder* for activity and transportation. Add 1ml sterile water to the tube, mix, and centrifuge a few seconds. Because multiple freeze-thaws may affect enzyme activity, store aliquots under -20°C.
4. Keep all of the reagents lids tightly capped when not in use to prevent evaporation, oxidation, and changes in pH.

## III. Principle:

The kit applies the unique buffer to rapidly lyse cell and inactivate cellular nuclease. DNA selectively adsorbs to silicified membrane in high salt solution. Cellular metabolite and proteins etc. are removed by serial of elution and centrifugation steps. Finally, low salt elution buffer elutes purified DNA from silica membrane.

## IV. Features:

1. Rapid and convenient. One preparation can be completed in 20 minutes.
2. The kit does not contain poisonous phenol and does not need a step of ethanol precipitation.

3. High-purity; unique membrane absorption and specialized washing for removing protein and other debris.
4. Multi-elution ensures high-purified DNA. The DNA yield is about 3-6  $\mu\text{g}$  from 200  $\mu\text{l}$  whole blood.

## V. Notes:

*Please read this section before your experiment.*

1. **All the centrifugation steps can be performed at room temperature.** Use a traditional centrifugal machine that the rotation speed can reach 13,000 rpm, such as Eppendorf 5415C.
2. The DNA yield is about 3-6  $\mu\text{g}$  from 200  $\mu\text{l}$  whole blood (the leukocyte number may vary in different samples, especially in patient blood, so the yield may vary largely).
3. Set water bath to 70°C before use.
4. For the best result, use fresh blood sample and avoid repeated freezing and thawing.

## VI. Procedure:

- *Add absolute ethanol to Buffer WB before use.*

1. Add 200  $\mu\text{l}$  fresh/cryogenic/anticoagulant and room temperature blood into 1.5 ml centrifuge tube.

*If the initial volume is less than 200  $\mu\text{l}$ , please add up to 200  $\mu\text{l}$  with Buffer BB; if the initial volume is between 200-300  $\mu\text{l}$ , increase the solution volume in the next step. If the initial volume is between 300  $\mu\text{l}$ —1 ml, please lyse erythrocyte first (as*

*described in appendix).*

2. Add 200 µl Buffer CB, shake for 15 seconds, and then add 20 µl protease K (20 mg/ml). Mix gently, and then incubate at 72°C water bath for 10 minutes. Solution should become clear.
3. Add 100 µl isopropanol, and then overturn to mix thoroughly. Flocculated precipitation may appear at this step.

*The proper strength and thorough mix is important for the DNA yield. Gentle vortex if necessary, but do not seriously agitate by hand to avoid shear DNA.*

4. Transfer the solution and flocculated precipitation into a Spin-column AC (Insert a spin-column AC into a collection tube), centrifuge at 10,000 rpm for 30 seconds, and discard the flow-through.
5. Add 500 µl Buffer IR and centrifuge at 12,000 rpm for 30 seconds. Discard the flow-through.
6. Add 700 µl Buffer WB **(please check if ethanol is added!)** and centrifuge at 12,000 rpm for 30 seconds. Discard the flow-through.
7. Add 500 µl Buffer WB and centrifuge at 12,000 rpm for 30 seconds. Discard the flow-through.
8. Put the Spin-column AC back to the collection tube and centrifuge at 13,000 rpm for 2 minutes. Remove rinsing buffer as much as possible. Ethanol will affect the down-stream reaction.
9. Transfer the Spin-column AC to a new microcentrifuge tube and add 100 µl preheated (65-70°C) Buffer EB. Place it at room temperature for 2-5 minutes. Centrifuge at

12,000 rpm for 1 minute. Add the flow-through onto the Spin-column AC and place it at room temperature for 2 minutes. Centrifuge at 12,000 rpm for 1 minute.

*The volume of elution buffer could be adjusted according to needs. Appropriate reduction of elution volume can increase concentration. However, the minimum volume is 20 µl. If the elution volume is less than 20 µl, elution efficiency and DNA yield can be affected.*

10. Store DNA at -20°C or apply to down-stream reactions.

## **VII. Appendix:**

(If sample volume is 300 µl-1 ml, lyse erythrocyte first)

1. Add 900 µl erythrocyte lysis buffer (not included; DP1801-A) to a 1.5 ml centrifuge tube or 3 ml erythrocyte lysis buffer to a 15 ml centrifuge tube.
2. Thoroughly mix the anticoagulant blood, add 300 µl/1 ml blood to 1.5 ml/15 ml centrifuge tube respectively, and invert for 6-8 times to make sure a thorough mix.
3. Place them at room temperature for 2-5 minutes.
4. Centrifuge at 12,000 rpm for 20 seconds (for 1.5 ml centrifuge tube)/ 2,000-3,000 rpm for 5 minutes (for 15 ml centrifuge tube) to remove the red supernatant, leave complete leukocyte mass and about 10 µl residual supernatant.
5. *White leukocyte mass (with a few erythrocytes) should be seen at the bottom of tube after centrifuge. If the most parts are red cell mass, it indicate that lysis of erythrocyte is not complete. Add 900 µl erythrocyte lysis buffer to a 1.5 ml centrifuge tube or 3 ml erythrocyte lysis buffer to a 15 ml centrifuge tube, resuspend the cell mass, and repeat*

*step 3-4 .*

6. Add 200 µl Buffer BB to suspend and fully disperse the leukocyte masses.
7. Isolate the blood genomic DNA as described in Procedures.

## VII. Troubleshooting:

| Problem                                                  | Possible reason                                                                                   | Advices                                                                                                                                |
|----------------------------------------------------------|---------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Blood clot appeared in sample                            | Improper storage of sample; Did not mix thoroughly or use improper anticoagulant collecting tube. | Discard the sample containing blood clot, re-collecting blood using anti-coagulating tube (containing EDTA, heparin, and citric acid). |
| Incomplete splitting of erythrocyte from white leukocyte | No adjusting to room temperature before sample splitting.                                         | Place sample to room temperature before use.                                                                                           |
|                                                          | Lysis time is not long enough.                                                                    | Prolonged to 15 minutes.                                                                                                               |
|                                                          | No thoroughly mix.                                                                                | Thoroughly mix the sample with the buffer.                                                                                             |
| Low DNA yield                                            | Low quantity of leukocyte in the sample.                                                          | Increase the initial quantity of blood.                                                                                                |
|                                                          | The storage time is too long.                                                                     | Use the fresh sample.                                                                                                                  |
|                                                          | Failure of protease K.                                                                            | The repeated melt and freeze of protease K should be avoided.                                                                          |
|                                                          | Not completely lysed; not thoroughly mix with isopropanol.                                        | Mix thoroughly after adding binding buffer and protease K. Mix thoroughly after add isopropanol, and then add into column.             |
|                                                          | Low elution efficiency.                                                                           | Make sure the correct operation of the steps 8 and 9; eluting only by EB.                                                              |
|                                                          |                                                                                                   |                                                                                                                                        |

|                                          |                                                                    |                                                                                                                                                                                                                                                                                                           |
|------------------------------------------|--------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                          | Some silicon based plasma membrane restrain endonuclease reaction. | Centrifuge genomic DNA solution at 13,000 rpm for 1 minute, carefully harvest supernatant.                                                                                                                                                                                                                |
| DNA size was less than 15 kb             | The blood sample is not fresh or improper stored.                  | Use fresh blood sample.                                                                                                                                                                                                                                                                                   |
|                                          | Incorrect operation leading to genome DNA shear.                   | Avoid vigorously mixing, transfer the samples by large diameter pipette tips.                                                                                                                                                                                                                             |
| A260/A280 too high                       | Silicon membrane interferes the value of A260/A280.                | Centrifuge genomic DNA solution at 13,000 rpm for 1 minute, carefully harvest supernatant.                                                                                                                                                                                                                |
| DNA remains slight colored after elution | Not enough washing.                                                | <ol style="list-style-type: none"> <li>1. After the step 7, if the flow-through is colored, repeat step 7 until the flow-through has no color.</li> <li>2. Use the eluted DNA as starting materials, and repeat the experiment again (skip the protease K digestion and 70°C incubation step).</li> </ol> |